

EXPERIMENT NO. 5

Quantitative Analysis

Read through the whole method before starting any practical work. Where appropriate, prepare a table for your results in the space provided.

Show your working and appropriate significant figures in the final answer to **each** step of your calculations.

- 1 In this experiment you will determine the percentage by mass of an impure sample of sodium hydrogencarbonate, NaHCO_3 . You will do this by titration with hydrochloric acid, HCl . The impurity in the sample is **X**. **X** is a sodium compound which does not react with HCl .

FB 1 is a mixture containing sodium hydrogencarbonate and **X**.

You are supplied with approximately 4 g of **FB 1**.

FB 2 is $0.105 \text{ mol dm}^{-3}$ hydrochloric acid, HCl .

methyl orange indicator

(a) Method

Preparing a solution of **FB 1**

- Weigh the 100 cm^3 beaker. Record the mass.
- Add between 2.8 g and 3.0 g of **FB 1** to the beaker.
- Reweigh the beaker with **FB 1**. Record the mass.
- Calculate and record the mass of **FB 1** used.
- Add approximately 50 cm^3 of distilled water to **FB 1** in the beaker.
- Stir the mixture with a glass rod until all the **FB 1** has dissolved.
- Transfer this solution into the 250 cm^3 volumetric flask.
- Wash the beaker with distilled water and transfer the washings to the volumetric flask.
- Add distilled water to the volumetric flask up to the mark.
- Shake the flask thoroughly.
- This solution of impure sodium hydrogencarbonate is **FB 3**. Label the flask **FB 3**.

Titration of **FB 3**

- Fill the burette with **FB 2**.
- Pipette 25.0 cm^3 of **FB 3** into a conical flask.
- Add approximately 3 drops of methyl orange indicator.
- Carry out a rough titration.
- Record your burette readings and the rough titre in the space below.

The rough titre is cm^3 .

- Carry out as many accurate titrations as you think necessary to obtain consistent results.
- Make sure any recorded results show the precision of your practical work.
- Record in a suitable form below all of your burette readings and the volume of **FB 2** added in each accurate titration.

I	
II	
III	
IV	
V	
VI	
VII	
VIII	

[8]

- (b) From your accurate titration results, obtain a suitable value for the volume of **FB 2** to be used in your calculations.
Show clearly how you obtained this value.

25.0 cm³ of **FB 3** required cm³ of **FB 2**. [1]

(c) Calculations

(i) Give your answers to (ii), (iii), (iv) and (v) to the appropriate number of significant figures. [1]

(ii) Calculate the number of moles of hydrochloric acid, HCl, in the volume of **FB 2** calculated in (b).

moles of HCl = mol [1]

(iii) Complete and balance the equation for the reaction of sodium hydrogencarbonate with hydrochloric acid. Include state symbols.



Deduce the number of moles of sodium hydrogencarbonate that reacted with the number of moles of HCl calculated in (ii).

moles of NaHCO₃ = mol [1]

(iv) Use your answer to (iii) to calculate the number of moles of sodium hydrogencarbonate in the **FB 1** that you weighed out.

moles of NaHCO₃ in **FB 1** used = mol [1]

(v) Calculate the percentage by mass of NaHCO₃ in **FB 1**. Ar of Na:23, H:1, C:12, O:16

percentage by mass of NaHCO₃ in **FB 1** = % [1]

[Total: 14]