## **EXPERIMENT NO. 3**

## 3 Qualitative Analysis

At each stage of any test you are to record details of the following.

- colour changes seen
- the formation of any precipitate
- the solubility of such precipitates in an excess of the reagent added

Where gases are released they should be identified by a test, **described in the appropriate place in your observations**.

You should indicate clearly at what stage in a test a change occurs. Marks are **not** given for chemical equations. **No additional tests for ions present should be attempted.** 

## If any solution is warmed, a boiling tube MUST be used.

Rinse and reuse test-tubes and boiling tubes where possible.

Where reagents are selected for use in a test, the name or correct formula of the element or compound must be given.

(a) **FB 4** and **FB 5** are solutions of salts each containing one cation and one anion from those listed in the Qualitative Analysis Notes. Carry out the following tests and record your observations in the table below.

	test	observations		
	lesi	FB 4	FB 5	
(i)	To a 1 cm depth of solution in a test-tube, add aqueous ammonia.	off-white ppt, turns brown on contact with air.	white ppt	
		ppt insoluble in excess	ppt insoluble in excess	
(ii)	To a 1 cm depth of solution in a test-tube, add a few drops of aqueous silver nitrate.	white ppt	no change	
(iii)	To a 1 cm depth of solution in a test-tube add a few drops of aqueous barium chloride or barium nitrate.	no Change	white ppt	

- (iv) Identify both ions in FB 4. cation  $M_n^{2+}$  anion  $C_{cay}$
- (v) Suggest the ions which may be present in **FB 5**.

cations  $M_{g^{2+}}/A_{caq}^{s+}$  anions  $SO_{4-}^{2-}/SO_{3-}^{2-}$ 

(vi) Select a reagent which could be used in a further test on **FB 5** to identify the **cation** present. Carry out your test and record your observations.

	test	observations				
	a 1 cm depth of <b>FB 5</b> in a st-tube, add	white ppt				
	ay. NaOH	PPt insoluble in excess				
The cation in <b>FB 5</b> is						

- (b) **FB 6** is a pale green salt containing two cations.
  - (i) What does this suggest about the identity of one of the cations in FB 6?

It contairs a transition metal ion.

Carry out the following tests and complete the table below.

		test	observations	
	(ii)	Place a spatula measure of <b>FB 6</b> in a hard-glass test-tube. Heat gently.	* solid dissolves to form a liquid * Vapours * solid turns brown	
	(iii)	Dissolve a small spatula measure of <b>FB 6</b> in a 2cm depth of distilled water in a test-tube. Use this solution for tests (iv) and (v).	solvid dissolves gives a yellow solution	
	(iv)	Pour about half the solution prepared in (iii) into a boiling tube and add aqueous sodium hydroxide, then	brawn/red-brawn ppt 19t insôlwble in excess	
		gently warm the mixture.	a colourless gas evolved which turns damp red litmus poper blue.	
	(v)	To the remainder of the solution prepared in (iii), add a few drops of aqueous potassium iodide, then	Solution darkens/ solution turns yellow-brown	
		add a few drops of starch solution.	turns blue-black	
(vi)	Identi	fy the cations present in <b>FB 6</b> .		
()			and	
(vii)	vada castin			
		JEWOK YEAU ON	[7]	

[Total: 14]