

## EXPERIMENT NO. 1

Metal carbonates react with dilute acids to produce carbon dioxide. You will identify the metal, **M**, in a metal carbonate,  $M_2CO_3$ , by measuring the volume of carbon dioxide produced during the reaction of  $M_2CO_3$  with excess hydrochloric acid.




**FA 2** is hydrochloric acid,  $HCl$ .

**FA 4** is  $M_2CO_3$ .

### (a) Method

Read **all** instructions before starting your practical work.

The diagrams below may help you in setting up your apparatus.



- Fill the tub with water to a depth of about 5 cm.
- Fill the 250 cm<sup>3</sup> measuring cylinder **completely** with water. Hold a piece of paper towel firmly over the top, invert the measuring cylinder and place it in the water in the tub.
- Remove the paper towel and clamp the inverted measuring cylinder so the open end is in the water just above the base of the tub.
- Use the 50 cm<sup>3</sup> measuring cylinder to place 50 cm<sup>3</sup> of **FA 2** into the reaction flask, labelled **X**.
- Check that the bung fits tightly in the neck of flask **X**, clamp flask **X**, and place the end of the delivery tube into the inverted 250 cm<sup>3</sup> measuring cylinder.
- Weigh the container with **FA 4** and record the mass in the space below.
- Remove the bung from the neck of the flask. Tip all the **FA 4** into the acid in the flask and replace the bung **immediately**. Remove the flask from the clamp and swirl it to mix the contents.
- Swirl the flask occasionally until no more gas is evolved. Replace the flask in the clamp.
- Reweigh the container and record the mass, and the mass of **FA 4** used, in the space below.
- When no more gas is collected, measure and record the final volume of gas in the measuring cylinder in the space below.

mass of tube + FA4 /g	21.70
mass of tube + residue /g	20.80
mass of FA4 used /g	0.90
initial volume of measuring cylinder /cm <sup>3</sup>	26
final volume of measuring cylinder /cm <sup>3</sup>	227
volume of CO <sub>2</sub> collected /cm <sup>3</sup>	201

[2]

**(b) Calculations**

Show your working and appropriate significant figures in the final answer to **each** step of your calculations.

- (i) Use the volume of gas you collected to calculate the number of moles of gas produced. [Assume that 1 mole of gas occupies 24.0 dm<sup>3</sup> under these conditions.]

$$\frac{201 \text{ cm}^3}{1000} = 0.201 \text{ dm}^3$$

$$\begin{array}{l} 24 \text{ dm}^3 \text{ --- } 1 \text{ mol} \\ 0.201 \text{ dm}^3 \text{ --- } x \end{array}$$

moles of gas =  $8.38 \times 10^{-3}$  mol

- (ii) Use your answer to (i) to deduce the number of moles of M<sub>2</sub>CO<sub>3</sub> used in the reaction.

mole ratio

$$\begin{array}{l} \text{M}_2\text{CO}_3 \quad : \quad \text{CO}_2 \\ 1 \quad \quad : \quad 1 \\ 8.38 \times 10^{-3} \quad : \quad 8.38 \times 10^{-3} \end{array}$$

moles of M<sub>2</sub>CO<sub>3</sub> =  $8.38 \times 10^{-3}$  mol

- (iii) Use your answer to (ii) and the mass of FA 4 used to calculate the relative formula mass, M<sub>r</sub>, of M<sub>2</sub>CO<sub>3</sub>.

$$\text{mol} = \frac{\text{mass}}{M_r}$$

$$8.38 \times 10^{-3} = \frac{0.90}{M_r} = 107.4$$

M<sub>r</sub> of M<sub>2</sub>CO<sub>3</sub> = 107.4

- (iv) Use your answer to (iii) and the Periodic Table to identify metal M. Explain your answer.

$$M_r \text{ of } \text{M}_2\text{CO}_3 = 107.4$$

$$M_r \text{ of } \text{CO}_3 = 60$$

$$M_r \text{ of } \text{M}_2 = 107.4 - 60 = 47.4$$

$$A_r \text{ of } \text{M} = \frac{47.4}{2} = \boxed{23.7}$$

M is Sodium - Na

As it is the nearest.

[4]

- (c) (i) A 250 cm<sup>3</sup> measuring cylinder can be read to ±1 cm<sup>3</sup>.

Calculate the maximum percentage error in your reading of the volume of gas.

$$\frac{1 \times 2}{201} \times 100$$

maximum percentage error = ...0.995... %

- (ii) It is likely that the volume of carbon dioxide that you collected was less than the theoretical volume.

Give **two** reasons why this volume is likely to be less than the theoretical volume.

In each case, suggest and explain a modification to the practical procedure that could help to reduce the difference in volume.

reason ... gas dissolves in water/solution ...

modification ... Use a gas syringe instead of measuring cylinder for gas collection ...

reason ... gas escapes before stopper inserted ...

modification ... Use more (excess) volume of a lower concentration of acid ...

[5]

[Total: 11]